Reference dimension in millimeters. Referencemått i millimeter. Referentiemåting in millimeter. Dimensions de référence en millimètres. Referenzabmessungen in Millimetern. Dimensioi indicate in millimeter. Referensial i millimeter. Dimensões de referência em milímetros. Dimensioes de referéncia em milímetros. Referensmått anges i millimeter.

1. **PRESERFLO™ MicroShunt GLAUCOMA DRAINAGE SYSTEM (English)**
2. **PRESERFLO™ MicroShunt GLAUKOMDRENSYSTEM (Danish)**
3. **PRESERFLO™ MicroShunt GLAUCOOM DRAINAGESYSTEEM (Dutch)**
4. **PRESERFLO™ MicroShunt SYSTÈME DE DRAINAGE POUR GLAUCOME (French)**
5. **PRESERFLO™ MicroShunt GLAUKOM-DRAINAGESYSTEM (German)**
6. **SISTEMA DI DRENAGGIO DEL GLAUCOMA PRESERFLO™ MicroShunt (Italian)**
7. **PRESERFLO™ MicroShunt DRENSYSTEM FOR GLAUKOM (Norwegian)**
8. **SISTEMA DE DRENAGEM PARA GLAUCOMA PRESERFLO™ MicroShunt (Portuguese)**
9. **SISTEMA DE DRENAJE PARA GLAUCOMA PRESERFLO™ MicroShunt (Spanish)**
10. **PRESERFLO™ MicroShunt DRÅNAGESYSTEM FÖR BEHANDLING AV GLAUKOM (Swedish)**

=====

PRESERFLO™ MicroShunt GLAUCOMA DRAINAGE SYSTEM (ENGLISH)
Instructions for Use

INTRODUCTION

The PRESERFLO™ MicroShunt (PFMS) is an implantable glaucoma drainage device made of an extremely flexible SIBS [poly(Syrene-block-Isobutylene-block-Styrene)] polymer with a tube of 350 µm outer diameter and a lumen of 70 µm. It has triangular fins that prevent migration of the tube into the AC. The device is designed to be implanted under the subconjunctival/Tenon space.

1. Surgeons should be well familiarized with the device with proper training.
2. Make sure to check the expiration date on each item.

CONTENTS

- PRESERFLO™ MicroShunt (1)	- 3mm Scleral Marker (1)	- 1mm Slit Angle Knife (1)
- Lasik Shields (3)	- Marker Pen (1)	- 25G Needle (1)

INDICATIONS FOR USE

The PRESERFLO™ MicroShunt Glaucoma Drainage System is intended for reduction of intraocular pressure in eyes of patients with primary open angle glaucoma where IOP remains uncontrollable while on maximum tolerated medical therapy and/ or where glaucoma progression warrants surgery.

CONTRAINDICATIONS

The implantation of the PRESERFLO™ MicroShunt is contraindicated if one or more of the following conditions exist: BACTERIAL CONJUNCTIVITIS; BACTERIAL CORNEAL ULCERS; ENDOPHTHALMITIS; ORBITAL CELLULITIS; BACTEREMIA OR SEPTICEMIA; ACTIVE SCLERITIS; UVEITIS; SEVERE DRY EYE; SEVERE BLEPHARITIS; PRE-EXISTING OCULAR OR SYSTEMIC PATHOLOGY THAT, IN THE OPINION OF THE SURGEON, IS LIKELY TO CAUSE POSTOPERATIVE COMPLICATIONS (E.G. SEVERE MYOPIA AND THIN CONJUNCTIVA) FOLLOWING IMPLANTATION OF THE DEVICE; PATIENTS DIAGNOSED WITH ANGLE CLOSURE GLAUCOMA.

WARNINGS

- Rx only. This device is restricted to sale by, or on the order of, a physician.
- For one-time use only. Do not reuse or re-sterilize. Reuse, or re-sterilization may compromise the structural integrity of the device and/ or lead to device failure which, in turn, may result in serious patient injury, illness, blindness or death. Reuse, or re-sterilization may also create a risk of contamination of the device and/ or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, blindness or death of the patient.
- After use, dispose of product and packaging in accordance with hospital, administrative, and/ or local government policy.
- Long term effects of Mitomycin C (MMC) with the use of this device have not been evaluated. Necessary precautions and interventions on the use of MMC are highly recommended.
- Viscoelasticity has not been tested with this device. However, in an emergency when all other therapies have failed, the use of hydroxyl methyl-cellulose (HPMC) may be an option. Use of HPMC should be a last resort to correct a flat chamber with the PRESERFLO™ MicroShunt and may risk loss of flow through the device for one or more weeks after use necessitating close or more frequent observation of IOP.
- Do not use on any patient who is too sensitive or allergic to metals.

PRECAUTIONS

- Should the PRESERFLO™ MicroShunt appear deformed, folded and/ or distorted, do not use.
- Avoid use of toothed forceps to handle device. McPherson type forceps are recommended.
- Do not use knife if it has been hit against or in contact with parts other than the eyeball. Such knife may be damaged and lose sharpness.
- The use of the PRESERFLO™ MicroShunt for secondary glaucoma has not been established.
- The device should not be subjected to direct contact with petrolatum-based (i.e. petroleum jelly) materials (e.g. ointments, dispersions, etc.).
- After-use discard with the blade stored inside the cover.

POTENTIAL COMPLICATIONS/ADVERSE EVENTS

The complications during and after surgery may include:

Glaucoma progression not controlled, difficulty in inserting the PFMS, extended surgical procedure, tube migration out of anterior chamber, flat anterior chamber, excessive bleeding in anterior chamber or eye, PFMS touches cornea or iris, intraocular pressure too high or low, viscoelastic used in anterior chamber, choroidal effusion or hemorrhage, retinal detachment, proliferative retinopathy, hypHEMA, hypotoni or hypotony maculopathy, phthisis bulbi, endophthalmitis, tube erosion through conjunctiva, tube block by iris or vitreous or fibrin, uveitis, diphlopia, aqueous misdirection, corneal complications (abrasion, edema, ulceration, infection, decompensation, bulbus keratopathy, endothelial cell loss, Descemet strip), partial or complete vision loss, globe perforation, bleb leak, blebitis, cystic bleb, bleb failure, pupillary block, ptosis, macular edema, prolonged inflammation, use of glaucoma medications, pain, conjunctival complications (dehiscence, dissection, hemorrhage, hypertonia, scar, ulcer), iris adhesions/synechia, cataract development or progression, explanation of the PFMS, encapsulation reaction, fibrin in anterior chamber, visual field damage, globe perforation, headache, vitreous hemorrhage, and suture related complications.

STERILIZATION

All PRESERFLO™ MicroShunt devices and accessories are sterile upon receipt. The PRESERFLO™ MicroShunt, scleral marker, and slit angle knife are sterilized by ETO sterilization per ISO 11135 guidelines and the marker pen, needle and Lasik Shields are sterilized by gamma radiation. Contents will remain sterile unless packaging is opened or damaged. Should the device or any accessory be inadvertently rendered non-sterile, or should the package be opened and/ or damaged, discard the device. **DO NOT RESTERILIZE.**

STORAGE

Store at controlled room temperature (20±5°C) and avoid excessive heat, excessive humidity, direct sunlight and water.

INSTRUCTIONS FOR USE

User must be properly trained on the use of the device prior to use.

Surgery:

1. **Anesthesia:** The type of anesthesia will be at the surgeon's discretion.

2. **Preparation of PRESERFLO™ MicroShunt:** Carefully examine the package containing the PRESERFLO™ MicroShunt for signs of damage that could compromise sterility. If damaged, discard the device. Remove PRESERFLO™ MicroShunt from sterile packaging onto sterile field. Wet the PRESERFLO™ MicroShunt using a solution of Balanced Salt Solution (BSS).
- Note: No cutting or modification of the PRESERFLO™ MicroShunt is allowed.**
3. **Preparation of Mitomycin C:** If the surgeon is applying MMC, follow the instructions for use for preparing 0.2-0.4mg/ml MMC. Prepare a syringe or line that can flush copious amounts of BSS (>2ml).
4. **Corneal Traction and Conjunctival Flap:** Perform corneal traction with 7-0 Vicryl suture at surgeon's discretion. In the planned surgical area, dissect a fimbria based subconjunctival/sub/Tenon pouch over a circumference of 90°-120° from near the limbus to at least 8mm deep.
5. **Application of MMC:** At the surgeon's discretion, apply 0.2-0.4mg/ml MMC soaked in 3 corneal lasik shield sponges under the subconjunctival pouch over a circumference of 90°-120° about 8mm deep for 2-3 minutes, remove the sponges and rinse with generous amounts of BSS (>2ml).

Note: Assure that the conjunctival wound edge does not come into contact with antimetabolite.

6. **Preparation of Anterior Chamber Entrance Site:** Assure there is no active bleeding in the area of the planned PFMS entry. Mark a point 3mm from the limbus utilizing 3mm scleral marker and marker pen. At the marked point, use a 1mm slit angle knife, create a shallow scleral flap at the mark just big enough to tuck the fin portion of the PRESERFLO™ MicroShunt. Make a 25G needle track into the anterior chamber (AC) by entering the scleral flap and advancing the needle so as to enter the angle just above the iris plane/MicroShunt.
7. **Insertion of PRESERFLO™ MicroShunt into Anterior Chamber:** Using forceps, position the the PRESERFLO™ MicroShunt tube with the bevel facing towards the cornea close to the needle track while keeping the fin parallel to the scleral surface. Gently insert the PRESERFLO™ MicroShunt into the track until the fin is tucked under the scleral flap. Use of excessive force will result in temporary bending of the PFMS and increase the difficulty of entrance. If not accessible, another tract can be made approximately 1mm to either side of the original tract. If the PRESERFLO™ MicroShunt is damaged during the procedure, use a new PRESERFLO™ MicroShunt and use a new tract.
- After successful insertion of the PRESERFLO™ MicroShunt, check that the PFMS tube position in the AC is not in contact with the cornea or iris and that there is flow of aqueous. If flow is not visible, use a 23g Thin wall 8mm bend cannula and inject BSS into the distal end of PRESERFLO™ MicroShunt to prime the tube. Confirm flow, then tuck the end of the PRESERFLO™ MicroShunt tube under the subconjunctival/sub/Tenon's flap on scleral surface making sure it is not caught in Tenon's. If there is no flow established, reposition or replace the device.
8. **Verify position of the MicroShunt in the anterior chamber angle:** It is recommended that intraoperative gonioscopy be performed to assess the position of the proximal tip of the PRESERFLO™ MicroShunt in the anterior chamber to ensure that there is sufficient clearance from the cornea or iris.

9. **Conjunctival Closure:** Reposition the Tenon's and conjunctiva to the limbus and suture using sutures with a long history of successful use by glaucoma specialists (i.e. nylon, vicryl). A moistened fluorescein strip is used to check for a conjunctival leak, which may be closed with additional sutures. Verify the presence of the proximal end of the PRESERFLO™ MicroShunt in the anterior chamber. Pull corneal traction suture out if it was used. Apply antibiotic and steroid medication post-op as appropriate.
10. After use, dispose of product and packaging in accordance with hospital, administrative, and/ or local government policy. Follow the MMC instructions for use for proper disposal of the MMC.

Post operatively:

1. Monitor the intraocular pressure at each follow-up to determine if the tube is open.
2. If the PRESERFLO™ MicroShunt is repositioned, removed and/ or replaced with another type of device based on the device performance, the conjunctiva should be cut at the limbus in a similar manner to the original procedure. The device should be exposed and the repositioning or removal conducted. If the device is removed, verify that there is no aqueous leakage from the track. If there is, use a Vicryl suture to seal the track. The limbus should then be re-sutured in a similar manner as the original surgery.

For any questions or concerns, contact: InnFocus Inc, USA; EMERGO EUROPE, Netherlands

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Docunet Corp.